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VICTOR CHONG, MD (OXFORD – UK)

1 TREATMENT PREPARATION: FOCUS


In order to guarantee a perfect laser beam focus, it is **important** to set the slit lamp ocular rings properly (diopter compensation adjustment). This **compulsory step** must be performed with the focusing rod of the slit lamp.

2 LASER SETTINGS

- Wavelength: **577 nm**
- Spot Size: **160 µm**
- SubLiminal® Mode / Duty Cycle: **5%**
- Exposure Time: **0.2 s**


Recommended Laser Lens: **Volk Area Centralis (0.94x)**

3 TITRATION PROCEDURE (COMPULSORY STEP BEFORE TREATMENT)

- A single spot is used to determine the thermal threshold of each patient.
- The power dose is evaluated in macular periphery in **a healthy area**.
- Increase the power level until a **barely visible threshold burn is observed***.
- Reduce the power to **50 % of the threshold power level** for treatment.

* While titrating, **do not exceed 1.2W power**. If no visible threshold happens at 1.2W, keep this power as a reference level and use half of it for treatment (i.e. 600mW).


4 USE OF THE MULTISPOT DELIVERY MODE TO IMPLEMENT THE SUBLIMINAL® TREATMENT

Macular Grid:

- Activation of the central fixation point.
- Resume® function activation.
- Set the spacing to 0 (laser impacts must be confluent).
- Deliver the burst of impacts in several times.

Squares: when the macular grid pattern cannot be used


- Resume® function activation.
- Set the spacing to 0 (laser impacts must be confluent).
- Deliver the burst of impacts in several times.

5 TREATMENTS

A. Treatment of Diabetic Macular Edema (DME)

OCT guided treatment: Treatment of identified edematous areas

Non-center involved diabetic macular edema:

The SubLiminal® treatment can be used alone. It induces biological changes in the RPE microenvironment that close the micro-aneurysms and dry the edema.

Foveal involved diabetic macular edema:


The best treatment option regarding this group of patients still remains an open question.

- In general, the thicker the retina, anti-VEGF therapy is more likely to be the first line treatment. Once the edema is settled, SubLiminal® laser can be used to reduce the number of injections.
- In patients with the fovea just slightly involved and/or without visual loss, SubLiminal® laser can be used as the first line treatment. If the edema deteriorates, then anti-VEGF therapy can be added.

B. Treatment of Central Serous Chorioretinopathy (CSCR):

ICGA / FA guided treatment

- ⇔ Treatment of the hyperfluorescent areas on mid-phase ICGA
- ⇔ Treatment of the "hot spots" on mid-phase FA

Important:

During the treatment:

- **No visible reaction** must be seen during the treatment.
- There is **no need to change power** with different degrees of edema.
- Laser impacts must be **confluent** (dense treatment) **and should not overlap**.
- **Transfoveal treatment is not recommended.**

After the treatment:

- DME Treatment Follow-up / **Results at 3 months.**
- CSCR Treatment Follow-up / **Results at 6 weeks.**

Refer to the **Supra Scan** user manual for detailed information on usage and safety measures.