1 INDICATIONS

- Acute angle-closure glaucoma
- Chronic angle-closure glaucoma
- Preventive treatment in eyes at risk for angle-closure glaucoma

Before performing an LPI, it is essential to perform a dynamic gonioscopy in order to differentiate between Appositional angle closure versus peripheral anterior synechiae. If the angle cannot be visualized by more than 200°, the LPI procedure is unlikely to be effective and could be a contraindication.

2 TREATMENT PREPARATION

>> FOCUS

In order to optimise laser beam focus, it is **important** to set the slit lamp oculars correctly for each user (diopter compensation adjustment).

This **compulsory step** must be performed using the focusing rod of the slit lamp.

>> LASER LENS

A contact lens such as an Abraham lens (+66D), Wise lens (+103D), Haag Streit CGIL or Volk Capsulotomy Lens should be used, under topical anesthesia. In order to ensure proper focus and minimize the risk of corneal epithelial burns, keep the eye open, limiting fine eye movements and apply an appropriate pressure in the case of bleeding.

3 TREATMENT SETTINGS

Clear iris:

- Offset: 0
- Energy: 2 - 3mJ

Dark iris:

- Offset: 0
- Energy: 3 - 5mJ

4 PRE-TREATMENT

Use Pilocarpine 2%, 30 minutes before the laser treatment, followed by apraclonidine 1% (this has the effect of iris vasoconstriction limiting the risk of bleeding and preventing IOP spikes).

Warning: Due to its mydriatic effect, apraclonidine must not be used without a pilocarpine pre-dose.

5 TREATMENT PROCEDURE

The YAG laser treatment can be performed **with or without** 532nm photocoagulation pre-treatment.


LPI can be done in one or two steps, the first is optional but preferable in the case of thick dark iris. Step one uses the 532nm laser photocoagulator to pre-treat the iris (around 5 impacts / power: 100-300mW; duration: 0,2-0,3s; spot size: 200µm).

The 532nm photocoagulation pre-treatment avoids the common issues associated with the use of pure YAG laser, including high energy levels and a high risk of iris haemorrhage.

It causes a thinning and contraction of the iris tissue (punched-out crater) which will be easier to puncture with low-energy YAG laser (combined technique).

Iridotomy site:

- The iridotomy site should be in the peripheral third of the iris.
- A crypt or a thinned area of the iris is recommended.
- 12 o'clock position is to be avoided (gas bubble).
- 11 o'clock and 1 o'clock positions are suitable (iridotomy covered by the lids avoiding diplopia), however the horizontal meridian position (slightly offset downwards) is favored as it avoids any interaction with the tear meniscus.

It is often possible to see aqueous humor flowing through the iridotomy site from the posterior chamber to the anterior chamber.

The optimal iridotomy size varies from 150 to 500 µm.

6 POST-TREATMENT

Apraclonidine 1% after the laser treatment to limit ocular hypertension.

Steroidal or non-steroidal anti-inflammatory drugs can be used with one drop of tropicamide.

Using a gonioscopy lens, please check for anterior peripheral synechial angle closures and / or plateau iris syndrome once the effects of the pilocarpine intake have ended.

The laser treatment may result in herpetic reactivation. A valaciclovir therapy is recommended for patients with history of ocular herpes.

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